



O-ring gantry mounted with a 6 MV, 500 cGy/minute linear accelerator on management studies and clinical trials involving the Exacta capable of rotating 360° around a patient during treatment, Verano system are eagerly awaited.

lowering total radiation dose exposure to pass-through normal tissues [4]. Treatment beams on the Vero are shaped by a multi-leaf collimator

comprised of 60 single-focused 5 mm leaves made of tungsten alloy capable of full over-center-travel. Vero has a minimum single treatment field size of 10x10 millimeters and a maximum treatment field size of

150x150 millimeters. What has the radiation oncology field buzz is that the linear accelerator and multi-leaf collimator are contained within a

Gimbal mount, allowing tracking of and treatment of mobile tumors. Further development of the Gimbal mount will expand the treatment

field size to 230x230 millimeters. Accurate tracking of tumors during quiet respiration by Vero ultimately may allow for lung treatment fields

to narrowly approximate PTV contours, adapt nearly in real-time with respiratory motion, and permit radiation dose intensification. Formal

References

1. Kunos C, DeBernardo R, Fabien J, Dobbins D, Zhang Y, et al. (2011) FDG-3 (7) & 7 GH¿QLWLRO RI FOLQLFDO WDUJHW YR
2. Kunos C (2013) Image-guided motion management. OMICS J Radiology 2: 410.
3. Lörich M, Rief N, Kühne M, Heckmann M, Fleckenstein J, et al. (2005) In vivo formation and repair of DNA double-strand breaks after computed tomography irradiations. Proc Natl Acad Sci U S A 102: 8984-8989.
4. Debuydt T, Penne R, Verellen D, Hrbacek J, Lang S, et al. (2012) Computer-aided motion management in radiotherapy. Phys Med Biol 57: 2997-3011.